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Correct application

- Ensure that only technically perfect, serviced and clean turbines, handpieces and contra-angles are used.
- The instruments must be inserted as far as possible.
- The instruments must be running before being placed on the surface.
- Avoid wedging or levering the instruments as this increases the risk of breakage.
- It is advisable to wear safety glasses, depending on the application.
- Excessive pressure must be avoided at all times.
- Bent or non-concentric rotary instruments and instruments with fractured and incorrectly shaped blades must be disposed of.

Cooling water spray

- To prevent undesirable heat build up during preparation, provide for sufficient cooling water spray (at least 50ml/min.).
- FG instruments with total length exceeding 22 mm or a head diameter of more than 2 mm require additional cooling.

Storage, disinfection, cleaning and sterilization

- Rotary instruments must be disinfected, cleaned and – whenever necessary – sterilized prior to first use on patients and immediately after use. A detailed description of the reprocessing of rotating instruments can be found in a separate chapter. The rotary instruments should be kept in its original packaging at room temperature and protected against dust and moisture until first use.
- Rotary instruments should be kept in hygienic stands, dishes or other suitable containers. The same applies to sterilized and sterile wrapped instruments. They must be protected against dust, moisture and recontamination during storage.
- Avoid contact with H₂O₂ (hydrogen peroxide). It attacks and damages tungsten carbide working parts, which curtails their service life.
- Avoid sterilization temperatures exceeding 180°C. Exceeding this temperature reduces the hardness of the working parts and curtails service life.
- Rotary instruments may corrode in a thermal disinfection unit. This may cause discolouration and curtail the service life.

Diamond instruments

For rotating diamond instruments with coarse diamond grain (green colour marking) and very coarse diamond grain (black colour marking), care should be taken for sufficient cooling and minimum contact pressure.

In order to achieve an optimal surface roughness it is essential to rework with a finisher after having used a coarse-grained or very coarse-grained diamond cutter.

Recommended speeds

- Maximum permissible speed have to be observed.

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Recommendations for reconditioning

For: Carbide Burs, Stainless Surgical Burs and Cutters, Stainless Root Canal Instruments and Diamond Instruments

In the case of the recommendations for reconditioning below, the manufacturer of the instruments has ensured that these are safe and suitable for use on the product. The user is responsible for ensuring that the reconditioning actually carried out with the equipment used, the materials and the personnel in the reconditioning facility achieves the desired aims and that the instruments are marked as required and that documentation is completed.

Rotary instruments must be disinfected, cleaned, dried and – wherever necessary - sterilised prior to first use on the patient and immediately after each use.

Limitations of reconditioning:

If the instruments are not identified as single use, all instruments can be used again. The end of serviceability is only determined by wear and tear or damage to the instruments. The dental surgeon is responsible for rejecting damaged and worn instruments at an early stage.

Storage:

Place used instruments immediately into a bath with special disinfection and cleaning solution. The method of use, reaction time and suitability of disinfectants for certain types of instruments are covered by the manufacturer's instructions, and it is essential that these are observed. We recommend only the use of preparations which have been tested and considered effective methods of disinfection by the German Institute for Hygiene and Microbiology (DGHM) or the Robert Koch Institute (RKI) (www.dghm.de and www.rki.de).

Reconditioning:

Reconditioning can be carried out by mechanically or manually. We recommend mechanical cleaning and disinfection; this procedure is validated. A standardised procedure is specified for manual reconditioning.

Mechanical cleaning:

The instruments should be rinsed under running water immediately before mechanical cleaning.

Place the instruments in suitable stands into a standard cleaning and disinfection machine.

Validated:

Cleaning and disinfection machine Miele G 7735 CD,

Vario TD programme:

- 4 minutes pre-washing with cold water
- 5 minutes washing at 55°C with an alkaline cleaning agent (0.5% Neodisher FA, Dr. Weigert, Hamburg)
- 3 minutes neutralisation with warm tap water (> 40°C)
- 2 minutes intermediate rinsing with warm tap water (>40°C)

Remove and dry instruments when programme is completed (preferably in accordance with the RKI recommendation using compressed air).

Carry out visual inspection for intactness and cleanliness, repeat the cleaning process if any contamination remains.

Manual cleaning:

Rinse the instruments under running water.

Heavily contaminated instruments should be cleaned first manually, by soaking the instruments in cold water and removing the particles of dirt with a soft brush.

Place the instruments in mesh baskets into an ultrasound bath with cleaning and disinfection liquid.

The method of use, reaction time and suitability of disinfectants for certain types of instruments are covered by the manufacturer's instructions, and it is essential that these are observed. We recommend only the use of preparations which have been tested and considered effective methods of disinfection by the German Institute for Hygiene and Microbiology (DGHM) or the Robert Koch Institute (RKI) as disinfection procedures found to be effective (www.dghm.de and www.rki.de).

At the end of the stipulated reaction time, rinse the instruments thoroughly with running water (preferably distilled water) and dry immediately (preferably in accordance with the RKI recommendation using compressed air).

Carry out visual inspection for intactness and cleanliness, repeat the cleaning process if any contamination remains. Since manual cleaning does not involve any automatic procedures, visual inspection is extremely important. The use of a magnifying device is recommended.

Sterilisation:

Suitable packaging should be selected for the instruments and the sterilisation procedure, and the instructions of the manufacturer of the sterilizer should be observed.

The validated procedure we recommend is steam sterilization with the following parameters:

- Fractioned pre-vacuum (3 times)
- Sterilization temperature 134°C
- Holding temperature: 5 minutes (full cycle)
- Drying time: 10 minutes

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Instructions for reconditioning non-stainless steel instruments

General directions:

Steel instruments (burs and cutters) made of non-stainless material cannot be sterilised in the autoclave. Steel instruments are therefore suitable only for applications for which the dental surgeon considers disinfection to be adequate.

In the case of the recommendations for reconditioning below, the manufacturer of the steel instruments has ensured that these are safe and suitable for use on the product. The user is responsible for ensuring that the reconditioning actually carried out with the equipment used, the materials and the personnel in the reconditioning facility achieves the desired aims and that the instruments are marked as required and that documentation is completed.

Steel instruments must be disinfected, cleaned and dried prior to first use on the patient and immediately after each use.

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Treatment after use:

Place used instruments immediately into a bath with special disinfection and cleaning solution. For steel instruments, a bath with anti-corrosion protection is required. We recommend the solution from DürrDental ID 220. The method of use, reaction time and suitability of disinfectants for certain types of instruments are covered by the manufacturer's instructions, and it is essential that these are observed. Only preparations should be used which have been tested and considered effective methods of disinfection by the German Institute for Hygiene and Microbiology (DGHM) or the Robert Koch Institute (RKI).

Reconditioning:

Reconditioning is manual in the case of steel instruments. Rinse the instruments thoroughly under running water.

Heavily contaminated instruments should be cleaned first with a brush. Place the instruments in mesh baskets into an ultrasound bath with cleaning and disinfection liquid.

For steel instruments, a bath with anti-corrosion protection is required. We recommend the solution from DürrDental ID 220. The method of use, reaction time and suitability of disinfectants for certain types of instruments are covered by the manufacturer's instructions, and it is essential that these are observed. Only preparations should be used which have been tested and considered effective methods of disinfection by the German Institute for Hygiene and Microbiology (DGHM) or the Robert Koch Institute (RKI) (www.dghm.de and www.rki.de).

At the end of the reaction time, instruments should immediately be dried thoroughly (preferably in accordance with the RKI recommendation using compressed air).

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